



Important Product Field Action Notification
Notice of Correction

Lyphochek Coagulation Control

Table with 10 columns: Level, Catalog Number, Lot Number, Lot Number, Lot Number, Lot Number, Lot Number, Lot Number, Lot Number, UDI. Rows include levels 1, 2, 3 and MiniPak with corresponding lot numbers and expiration dates.

This Notice is intended for the end-users of this product.
If you are not the end-user, please forward this Notice to the appropriate laboratory personnel.

Date: January 31, 2023
To: Clinical Laboratory Director, Manager or Supervisor
From: Bio-Rad Laboratories QSD Regulatory Affairs
Product: Lyphochek Coagulation Control
Issue: We have received intermittent complaints about flakes/clots in Lyphochek Coagulation Control. This product is prepared from non-defibrinated plasma; sporadic flakes/clots can form when the reconstituted product is frozen. We are, therefore, removing the reconstituted frozen and freeze-thaw cycles claims and are revising the LIMITATIONS section to instruct customers to discard the vials if they find evidence of fibrin clots.

The STORAGE AND STABILITY and LIMITATIONS sections in the Instructions for Use (IFU) have been revised as shown below. (This is a permanent change and will also be reflected in the IFU of all future lots).

All other claims for Lyphochek Coagulation Control remain unchanged.

Revised Claims



Significant changes are highlighted!

STORAGE AND STABILITY

After reconstituting and storing tightly capped at 2 to 25°C, this product will be stable as follows:

- All analytes: 48 hours

LIMITATIONS

If there is evidence of microbial contamination, excessive turbidity, or fibrin clots in the reconstituted product, discard the vial.

Actions Required

- Please obtain the revised (Rev. 1) IFU through the Internet at http://myeinserts.qcnet.com. Discard all previous versions that you may have at your location.
Contact Bio-Rad Laboratories Technical Support Department (USA) at 1-800-854-6737, option #2, or by email at qsd.techservice@bio-rad.com if you have technical questions regarding this notice.

We appreciate your patience and apologize for any inconvenience this may have caused.