

Carestream Health Inc.
150 Verona Street
Rochester, NY, 14608

Date: September 20, 2023

MEDICAL DEVICE FIELD CORRECTIVE ACTION

To: Director of Radiology/Diagnostic Imaging, Radiology Administrators, and Radiology/Diagnostic Imaging Managers

You are receiving this communication to inform you that Carestream is conducting a Field Corrective Action involving Carestream DRX-Revolution Mobile X-ray System(s) in your possession.

Issue Description

There is a potential for unexpected failure of the electrical components within the Carestream Health Inc. generator installed in your DRX-Revolution Mobile X-ray System(s). This can lead to temporary and self-contained thermal overload within the generator. Should an event like this occur, the system will become inoperable and loud noise, burnt smell and smoke may be detectable.

The historical rate of incidence is small and predictive of the rate of future incidents. Carestream has determined that such an event is unlikely to occur.

Actions To Be Taken By User Facility:

- Inform all personnel that utilize the system of the potential issue
- Should an event of this nature occur:
 - Move the system away from patients/staff.
 - The system will be inoperable.
 - Call the Carestream Customer Care Center in the U.S. at 1-800-328-2910, your local Carestream Service support number or Carestream Health Authorized Representative to place a service call.
- Complete the "Consignee Notification Acknowledgement" form and respond as instructed on the form.

Corrective Action To Be Taken By Carestream:

You will be contacted by a Carestream Service Engineer or Carestream Health Authorized Representative to schedule a time for your system(s) to be serviced. At the scheduled time, a Carestream Service Engineer or Carestream Health Authorized Representative will inspect your DRX-Revolution Mobile X-ray System and replace generator component(s) which will resolve the identified Carestream generator issue.


If you have any questions or concerns please contact the Carestream Customer Care Center in the U.S. at 1-800-328-2910, available 7 days per week on a 24-hour basis. Outside of the U.S., please call your local Carestream Service support number.

If you have distributed the device outside your facility, please alert your customer(s) of this Field Corrective Action and contact the Carestream Customer Care Center as instructed above.

This Field Corrective Action is being made with the knowledge of the US Food and Drug Administration.

We regret any inconvenience this may cause to your operations.

Kind Regards,



Glen Nicholson
Product Line Manager
Mobile Solutions
Carestream Health Inc.