

Attachment A: Serial No of Medical Device

Batch No.
RLA071696G
RLA148346G
RLA148317G
RLA148278G
RLA142203G
RLA142236G
RLA130029G
RLA150247G
RLA150244G
RLA129364G
RLA150252G
RLA150784G
RLA150237G
RLA150229G
RLA151372G
RLA112744G
RLA143013G
RLA150227G
RLA150246G
RLA150217G
RLA150224G
RLA130058G
RLA143011G
RLA148276G
RLA124304G
RLA130030G
RLA148332G
RLA143050G
RLA114510G
RLA094723G
RLA129355G
RLA130051G
RLA142227G
RLA081166G
RLA119033G
RLA130041G
RLA130049G
RLA119022G
RLA089111G
RLA129361G
RLA090305G
RLA102131G
RLA090314G
RLA090657G
RLA112814G
RLA107231G
RLA114518G
RLA088778G
RLA088850G

RLA114853G
RLA080193G
RLA112743G
RLA112736G
RLA112738G
RLA102183G
RLA107240G
RLA105419G
RLA094722G
RLA102545G
RLA107241G
RLA094716G
RLA094721G
RLA107238G
RLA102542G
RLA094720G
RLA102539G
RLA102187G
RLA102130G
RLA081164G
RLA094718G
RLA094719G
RLA094724G
RLA060024G
RLA094166G
RLA061152G
RLA094156G
RLA081363G
RLA082513G
RLA081170G
RLA080194G
RLA081714G
RLA082511G
RLA059804G
RLA061154G
RLA081365G
RLA071015G
RLA060078G
RLA072446G
RLA061182G
RLA080192G
RLA060068G
RLA071018G
RLA071017G
RLA009773G
RLA071976G
RLA059805G
RLA061177G
RLA031190G
RLA059815G

RLA052181G
RLA034751G
RLA048261G
RLA052183G
RLA052184G
RLA048340G
RLA034763G
RLA047541G
RLA046516G
RLA048263G
RLA051894G
RLA048257G
RLA013021G
RLA034728G
RLA034755G
RLA034760G
RLA047535G
RLA034735G
RLA032981G
RLA034726G
RLA028604G
RLA034754G
RLA025246G
RLA027804G
RLA034729G
RLA013018G
RLA027792G
RLA021771G
RLA009754G
RLA031106G
RLA023537G
RLA025797G
RLA028530G
RLA017517G
RLA023404G
RLA017516G
RLA023401G
RLA023400G
RLA023403G
RLA015918G
RLA015913G
RLA017332G
RLA015932G
RLA017512G
RLA014921G
RLA014998G
RLA015031G
RLA015927G
RLA014926G
RLA013040G
RLA009628G

RLA215327G
RLA262520G
RLA224231G
RLA249877G
RLA222815G
RLA238171G
RLA234679G
RLA215336G
RLA222912G
RLA222909G
RLA249071G
RLA222808G
RLA224230G
RLA208841G
RLA215330G
RLA215341G
RLA213065G
RLA215329G
RLA208774G
RLA213068G
RLA213130G
RLA208842G
RLA215328G
RLA215332G
RLA208496G
RLA207772G
RLA207779G
RLA207780G
RLA207789G
RLA195928G
RLA208495G
RLA207781G
RLA195930G
RLA195629G
RLA195931G
RLA180826G
RLA190065G
RLA191932G
RLA195934G
RLA207705G
RLA180821G
RLA187496G
RLA187513G
RLA190062G
RLA180822G
RLA187483G
RLA178435G
RLA164300G
RLA178440G
RLA164299G
RLA159803G

APPENDIX B: Customer Communication letter



Medtronic International, Ltd. (Singapore Branch)

(Co. Reg. No. S98FC5604C)
50 Pasir Panjang Road
#04-51
Mapletree Business City
Singapore 117384
www.medtronic.com

tel 65.6870.5300
fax 65.6482.0300

Medical Device Correction

**Reveal LINQ™ with TruRhythm™ Insertable Cardiac Monitoring Systems
Brady & Pause Detections Disabled Following Partial Electrical Reset**

2nd June 2021

Attention: Risk management Director and O.R Materials Management

CC: The Chairman Medical Board and relevant Head of Departments

Affected Model Number	Model Description
LNQ11	Reveal LINQ with TruRhythm

2nd June 2021

Dear Risk Manager or Healthcare Professional,

Please share this notification with the Cardiology and cardiac monitoring departments, Pacemaker/Device clinic leadership, and physicians who implant or manage patients with Reveal LINQ™ with TruRhythm™ insertable cardiac monitors (ICMs).

This letter is to inform you that Reveal LINQ with TruRhythm ICMs that undergo a partial electrical reset appear to be programmed "ON," but are no longer able to detect and report Brady and Pause events to clinicians. Medtronic estimates that **0.049%** of Reveal LINQ with TruRhythm ICMs have experienced a partial electrical reset resulting in the inability to detect Brady and Pause events. While there is a potential for underreporting due to lack of awareness that an electrical reset has occurred, there have been zero (0) serious or permanent harms or deaths reported as a result of this issue. After a partial electrical reset, these Brady and Pause episode types will not be reported to the clinician.

- Currently implanted/distributed Reveal LINQ with TruRhythm ICMs will receive a future software update to correct this issue delivered via the Model 2090 and Encore™ programmers. The corrective fix is anticipated to be available. Availability of the software will be communicated once Medtronic has obtained the necessary regulatory approvals.

- There will be an update for future manufactured Reveal LINQ with TruRhythm ICMs, which is anticipated to be available beginning late calendar year 2021. Medtronic will inform physicians once this manufacturing update is implemented into newly manufactured Reveal LINQ with TruRhythm ICMs.

ISSUE DESCRIPTION

Medtronic has identified that Reveal LINQ with TruRhythm ICMs that undergo a partial electrical reset appear to be programmed "ON," but are no longer able to detect and report Brady and Pause events. A partial electrical reset is normal behavior that can occur when the device detects a possible issue with the device software. However, an error in the partial electrical reset implementation is causing this unintended behavior.

All Reveal LINQ with TruRhythm ICMs currently in distribution are susceptible to this issue. Through 10 May 2021, Medtronic has received 87 complaints related to an electrical reset. The projected rate of a Reveal LINQ with TruRhythm ICM experiencing a partial electrical reset that results in the inability to detect Brady and Pause events is 0.056% at 36 months. Complaint data suggests the majority of electrical resets were associated with Electromagnetic Interference (EMI) due to cardioversion or electrocautery. Potential harms include those associated with the risk of a delayed medical intervention or missed diagnosis for Brady and Pause events, and an explant procedure.

If a partial electrical reset occurs, CareLink™, Model 2090 and Encore programmer software and Reveal LINQ™ Mobile Manager (LMM) will continue to indicate that detection parameters are "ON;" however, Brady and Pause events will not be automatically collected. The Patient Assistant (Patient Activator) will continue to function to manually trigger ECG collection, store the tracing, and mark symptoms.

Tachy and AT/AF detections are not affected by a partial electrical reset.

HOSPITAL RISK MANAGER ACTIONS

1. Please share this notification with the Cardiology and cardiac monitoring departments, Pacemaker/Device Clinic leadership, and physicians who implant or manage patients with LINQ II insertable cardiac monitors (ICMs).
2. Complete the enclosed Confirmation Form and provide to your Medtronic representative.

PATIENT MANAGEMENT RECOMMENDATIONS

All patients, including those on CareLink, should be carefully monitored for reports of an electrical reset condition. Follow instructions below.

- During in person or remote follow-up: If a device experiences an electrical reset, clinicians will be informed via programmer pop-up or CareLink display message. Actively monitor for these notifications at each patient follow-up, and contact Medtronic Representative should you receive an alert.
Note: Once cleared, electrical reset notifications are no longer accessible.
- Retroactively: Review the Brady lifetime episode counter from the most recent session report (CareLink or in-office). If a report is not available, consider scheduling a follow-up for each patient being monitored for Brady or Pause events. Review the Brady lifetime episode counter:
 - If the lifetime count for Brady is non-zero, a partial electrical reset has not occurred.
 - If the lifetime count for Brady is zero, and the Brady detection parameter indicates it is "ON," a partial electrical reset may have occurred. Contact your Medtronic Representative.

Patients with a confirmed partial electrical reset:

- Medtronic medical staff, in consultation with our Independent Physician Quality Panel, recommends against device replacement for patients being monitoring for Tachy or AT/AF; continue normal patient follow-up.
- When monitoring for Brady or Pause events, it is important to note that the Patient Assistant (Patient Activator) will continue to manually mark symptoms even after a partial electrical reset. Patient-activated recordings are not impacted by this issue. If patients require monitoring for Brady and/or Pause events, and it is not acceptable to wait for the software update to become available (see details below), consider device replacement. Recognize that exposure to EMI could introduce this issue for new device implants that occur before the manufacturing update is implemented anticipated late calendar year 2021.
- As a reminder, per the Reveal LINQ with TruRhythm ICM's Instructions for Use, contact your Medtronic representative anytime an electrical reset occurs.

FUTURE SOFTWARE UPDATE AVAILABILITY

Medtronic is developing a programmer-delivered software update to correct this issue for Reveal LINQ with TruRhythm ICMs currently implanted or in distribution. Anticipated availability is early calendar year 2022 ; Medtronic representatives will inform you of the availability and work with you to install the software onto clinic and hospital 2090 and Encore programmers. LMM application software will be unable to deliver the software update for this issue. In order for patients with Reveal LINQ with TruRhythm ICMs to receive the update, the device will need to be interrogated with an updated 2090 or Encore programmer.

Medtronic will notify all applicable regulatory agencies about this matter. Please share this notification with others in your organization as appropriate.

Please notify your Medtronic representative of any adverse events or quality problems associated with your use of this product.

We sincerely regret any difficulties this may cause you and your patients. Medtronic remains dedicated to patient safety and will continue to monitor device performance to ensure we meet your needs and those of your patients

Sincerely,



Sincerely,
Diana Teo
Medtronic QRA Head
Singapore and Malaysia



Sincerely,
Chloe Tan
Medtronic QRA Head
Indochina and Frontier Market
Plus



Sincerely,
Parichart Bunjobchokchai
Medtronic QRA Senior Manager
Thailand

Enclosure: Customer Confirmation Form



Medtronic International, Ltd. (Singapore Branch)

(Co. Reg. No. S98FC5604C)
50 Pasir Panjang Road
#04-51
Mapletree Business City
Singapore 117384
www.medtronic.com

tel 65.6870.5300
fax 65.6482.0300

Customer Confirmation Form

ALL CUSTOMERS PLEASE COMPLETE THE FORM IN ITS ENTIRETY

Customer Contact Details	Medtronic Contact Details
Physician / HCP/Distributor :	Name:
Address:	Contact:
Phone no:	Email:
E-mail:	

For completion by Medtronic Customers Only – Please complete all fields below and return immediately

In the event you no longer implant and/or manage patients with Reveal LINQ™ please provide a detailed explanation in the space below so that Medtronic’s records can be updated accordingly. Thank you

Note: The addressee may continue to receive reminders of this notice until a response is received.

By signing this form, I confirm that I have read the Medical Device Correction Notification Letter, listed below dated 2nd June 2021 from Medtronic and taken appropriate action for Reveal LINQ™ with TruRhythm™ Insertable Cardiac Monitoring Systems - Brady & Pause Detections Disabled Following Partial Electrical Reset

Name: _____ (print) Signature: _____ Stamp: _____ Date: _____