



**PRODUCT NOTIFICATION**

**IDS-23-4851-FA**

**BD BBL™ Sensi-Disc™**

06 February 2024

Dear Customer,

CC: Chairman Medical Board and relevant head of department

**Type of Field Action:**

- Removal (Recall): For Catalog #291270 only (Impacted Lot: 3093584 & 3184136)
- Product Advisory (No action required): All other Catalogs listed in Appendix 1

**Affected Product:**

Refer to Appendix 1: List of Affected Products

**Description of the Problem:**

BD confirmed through a recent evaluation of BD BBL™ Sensi-Disc™ product that twenty-eight (28) out of thirty (30) antimicrobial discs showed reproducibility, accuracy, and/or Quality Control (QC) failures when tested with *Haemophilus* spp. Testing with other bacterial pathogens (where indicated) is not affected.

**Clinical Risk Statement:**

Based on findings from internal and reference lab testing, there is a possibility of reproducibility, accuracy and/or QC failures in antibiotic susceptibility testing (AST) for *H. influenzae*. Performance is highly variable depending on plate manufacturer, AST guidelines utilized, and the antibiotic tested. This may cause product discard, delayed results, or additional adverse diagnostic outcomes, such as a delay in diagnosis, the selection of inappropriate antibiotics or extended duration of antibiotic exposure and the treatment process.

**Complaint & Adverse Event Statement:**

To date, there have been no adverse events worldwide related to this issue.

**Actions for Clinical Users:**

1. Users should refrain from conducting AST for *H. influenzae* with impacted Sensi-Discs and cease use of catalog number 291270, which is solely intended for AST for *H. influenzae*, to mitigate these possible risks. The affected product can continue to be used with other bacterial pathogens (where indicated).
2. It is not necessary to review previous test results and no additional clinical actions are recommended.

**Action to be taken by BD:**

1. BD is investigating root cause and will implement appropriate corrective and preventative measures to prevent recurrence.
2. BD will issue credit to customers that have catalog number 291270, following receipt of the completed Customer Response Form.



3. Future lots of BD BBL™ Sensi-Disc™ provided (excluding catalog number 291270) will be labeled with a sticker notifying customers “This product should not be used for the semi-quantitative in vitro susceptibility testing of *Haemophilus influenzae*”.

**Please Take the Following Actions:**

1. Immediately inspect your inventory for the specific catalog and lot numbers listed to determine if the BD BBL™ Sensi-Disc™ in your possession are impacted. Ensure the contents of this notification are read and understood.
2. If you have catalog number 291270 in your inventory, return it to your distributor for further instructions.
3. Other catalog numbers with lot numbers in scope of this action can continue to be used with other bacterial pathogens (where indicated).
4. Share and post this customer letter within your facility network and forward to any customers to whom you may have distributed the product to ensure awareness.
5. Complete the attached Customer Response Form and return it to the BD contact noted on the form indicating whether or not you have any of the affected lots so that BD may acknowledge your receipt of this notification.

**Appendix:**

Appendix 1: List of Affected Products.

BD is committed to advancing the world of health. Our primary objectives are patient and user safety and providing you with quality products. We apologize for any inconvenience this issue may have caused you and thank you in advance for helping us to resolve this matter as quickly and effectively as possible.

Yours Sincerely,

DocuSigned by:  
*Karena Han*  
Signer Name: Karena Han  
Signing Reason: I approve this document  
Signing Time: 06-Feb-2024 | 3:35:10 AM GMT  
992DAF448F77410697FB3524C491D633

---

Karena Han  
Quality Manager, Southeast Asia



**Appendix 1**  
**List of Affected Products**

Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
230705	Sensi Disc Ampicillin - 10 µg	1334111	(01) 0 038290 230705 1	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
230705	Sensi Disc Ampicillin - 10 µg	2025840	(01) 0 038290 230705 1	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
230705	Sensi Disc Ampicillin - 10 µg	2146318	(01) 0 038290 230705 1	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
230705	Sensi Disc Ampicillin - 10 µg	2242344	(01) 0 038290 230705 1	30-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
230705	Sensi Disc Ampicillin - 10 µg	2276051	(01) 0 038290 230705 1	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
230705	Sensi Disc Ampicillin - 10 µg	2333100	(01) 0 038290 230705 1	31-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
230998	Sensi Disc Tetracycline – 30 µg	1242880	(01) 0 038290 230998 7	30-Sep-2026	single (one cartridge) or Ten Pack (10 cartridges)
230998	Sensi Disc Tetracycline – 30 µg	1277918	(01) 0 038290 230998 7	31-Oct-2026	single (one cartridge) or Ten Pack (10 cartridges)
230998	Sensi Disc Tetracycline – 30 µg	1301570	(01) 0 038290 230998 7	30-Nov-2026	single (one cartridge) or Ten Pack (10 cartridges)
230998	Sensi Disc Tetracycline – 30 µg	2025088	(01) 0 038290 230998 7	28-Feb-2027	single (one cartridge) or Ten Pack (10 cartridges)
230998	Sensi Disc Tetracycline – 30 µg	2056233	(01) 0 038290 230998 7	31-Mar-2027	single (one cartridge) or Ten Pack (10 cartridges)
230998	Sensi Disc Tetracycline – 30 µg	2118569	(01) 0 038290 230998 7	31-May-2027	single (one cartridge) or Ten Pack (10 cartridges)
230998	Sensi Disc Tetracycline – 30 µg	2276681	(01) 0 038290 230998 7	31-Oct-2027	single (one cartridge) or Ten Pack (10 cartridges)
230998	Sensi Disc Tetracycline – 30 µg	3094187	(01) 0 038290 230998 7	30-Apr-2028	single (one cartridge) or Ten Pack (10 cartridges)
230998	Sensi Disc Tetracycline – 30 µg	9009914	(01) 0 038290 230998 7	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
230998	Sensi Disc Tetracycline – 30 µg	9046895	(01) 0 038290 230998 7	31-Mar-2024	single (one cartridge) or Ten Pack (10 cartridges)
230998	Sensi Disc Tetracycline – 30 µg	0065366	(01) 0 038290 230998 7	31-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
230998	Sensi Disc Tetracycline – 30 µg	0178210	(01) 0 038290 230998 7	31-Jul-2025	single (one cartridge) or Ten Pack (10 cartridges)
231250	Sensi Disc Minocycline – 30 µg	1085656	(01) 0 038290 231250 5	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231250	Sensi Disc Minocycline – 30 µg	1277053	(01) 0 038290 231250 5	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
231250	Sensi Disc Minocycline – 30 µg	1301694	(01) 0 038290 231250 5	30-Nov-2024	single (one cartridge) or Ten Pack (10 cartridges)
231250	Sensi Disc Minocycline – 30 µg	1334451	(01) 0 038290 231250 5	31-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
231250	Sensi Disc Minocycline – 30 µg	2025108	(01) 0 038290 231250 5	28-Feb-2025	single (one cartridge) or Ten Pack (10 cartridges)
231250	Sensi Disc Minocycline – 30 µg	2119580	(01) 0 038290 231250 5	31-May-2025	single (one cartridge) or Ten Pack (10 cartridges)



2 International Business Park Road  
 The Strategy #08-08  
 Singapore 609930  
 Registration No. 201114149N  
**bd.com**

Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
231250	Sensi Disc Minocycline – 30 µg	2146396	(01) 0 038290 231250 5	30-Jun-2025	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	1060689	(01) 0 038290 231536 0	31-Mar-2024	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	1116935	(01) 0 038290 231536 0	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	1210905	(01) 0 038290 231536 0	31-Aug-2024	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	1242939	(01) 0 038290 231536 0	30-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	1306103	(01) 0 038290 231536 0	30-Nov-2024	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	1334461	(01) 0 038290 231536 0	31-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	2004540	(01) 0 038290 231536 0	31-Jan-2025	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	2025116	(01) 0 038290 231536 0	28-Feb-2025	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	2056456	(01) 0 038290 231536 0	31-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	2108378	(01) 0 038290 231536 0	30-Apr-2025	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	2146453	(01) 0 038290 231536 0	30-Jun-2025	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	2301307	(01) 0 038290 231536 0	30-Nov-2025	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	3003869	(01) 0 038290 231536 0	31-Jan-2026	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	3062089	(01) 0 038290 231536 0	31-Mar-2026	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	3094151	(01) 0 038290 231536 0	30-Apr-2026	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	3118079	(01) 0 038290 231536 0	31-May-2026	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	3152968	(01) 0 038290 231536 0	08-Jun-2026	single (one cartridge) or Ten Pack (10 cartridges)
231536	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	0343394	(01) 0 038290 231536 0	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231541	Sensi Disc Rifampin – 5 µg	1278110	(01) 0 038290 231541 4	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
231541	Sensi Disc Rifampin – 5 µg	1334479	(01) 0 038290 231541 4	31-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
231541	Sensi Disc Rifampin – 5 µg	2004554	(01) 0 038290 231541 4	31-Jan-2025	single (one cartridge) or Ten Pack (10 cartridges)
231541	Sensi Disc Rifampin – 5 µg	2025125	(01) 0 038290 231541 4	28-Feb-2025	single (one cartridge) or Ten Pack (10 cartridges)
231541	Sensi Disc Rifampin – 5 µg	2056466	(01) 0 038290 231541 4	31-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
231541	Sensi Disc Rifampin – 5 µg	2242510	(01) 0 038290 231541 4	30-Sep-2025	single (one cartridge) or Ten Pack (10 cartridges)



Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
231541	Sensi Disc Rifampin – 5 µg	3214931	(01) 0 038290 231541 4	01-Aug-2026	single (one cartridge) or Ten Pack (10 cartridges)
231606	Sensi Disc Cefotaxime – 30 µg	2004580	(01) 0 038290 231606 0	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
231606	Sensi Disc Cefotaxime – 30 µg	2025164	(01) 0 038290 231606 0	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231606	Sensi Disc Cefotaxime – 30 µg	2090467	(01) 0 038290 231606 0	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231606	Sensi Disc Cefotaxime – 30 µg	2146516	(01) 0 038290 231606 0	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
231606	Sensi Disc Cefotaxime – 30 µg	2180432	(01) 0 038290 231606 0	31-Jul-2024	single (one cartridge) or Ten Pack (10 cartridges)
231606	Sensi Disc Cefotaxime – 30 µg	2301360	(01) 0 038290 231606 0	30-Nov-2024	single (one cartridge) or Ten Pack (10 cartridges)
231606	Sensi Disc Cefotaxime – 30 µg	3003886	(01) 0 038290 231606 0	31-Jan-2025	single (one cartridge) or Ten Pack (10 cartridges)
231628	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	2146520	(01) 0 038290 231628 2	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231628	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	2210201	(01) 0 038290 231628 2	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231628	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	2276143	(01) 0 038290 231628 2	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231628	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	3062095	(01) 0 038290 231628 2	30-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
231634	Sensi Disc Ceftriaxone - 30 µg	1033870	(01) 0 038290 231634 3	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231634	Sensi Disc Ceftriaxone - 30 µg	1244758	(01) 0 038290 231634 3	30-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
231634	Sensi Disc Ceftriaxone - 30 µg	2026381	(01) 0 038290 231634 3	28-Feb-2025	single (one cartridge) or Ten Pack (10 cartridges)
231634	Sensi Disc Ceftriaxone - 30 µg	2056573	(01) 0 038290 231634 3	31-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
231634	Sensi Disc Ceftriaxone - 30 µg	2277488	(01) 0 038290 231634 3	31-Oct-2025	single (one cartridge) or Ten Pack (10 cartridges)
231644	Sensi Disc Imipenem – 10 µg	1335006	(01) 0 038290 231644 2	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231644	Sensi Disc Imipenem – 10 µg	2056580	(01) 0 038290 231644 2	31-Mar-2024	single (one cartridge) or Ten Pack (10 cartridges)
231644	Sensi Disc Imipenem – 10 µg	2091918	(01) 0 038290 231644 2	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231644	Sensi Disc Imipenem – 10 µg	2210318	(01) 0 038290 231644 2	31-Aug-2024	single (one cartridge) or Ten Pack (10 cartridges)
231644	Sensi Disc Imipenem – 10 µg	2301385	(01) 0 038290 231644 2	30-Nov-2024	single (one cartridge) or Ten Pack (10 cartridges)
231644	Sensi Disc Imipenem – 10 µg	3062101	(01) 0 038290 231644 2	31-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
231644	Sensi Disc Imipenem – 10 µg	3094134	(01) 0 038290 231644 2	30-Apr-2025	single (one cartridge) or Ten Pack (10 cartridges)
231673	Sensi Disc Cefpodoxime – 10 µg	3109152	(01) 00382902316732	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)



2 International Business Park Road  
 The Strategy #08-08  
 Singapore 609930  
 Registration No. 201114149N  
[bd.com](http://bd.com)

Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
231673	Sensi Disc Cefpodoxime – 10 µg	3157106	(01) 00382902316732	08-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
231673	Sensi Disc Cefpodoxime – 10 µg	3241174	(01) 00382902316732	27-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
231691	Sensi Disc Piperacillin/Tazobactam – 100/10 µg	2026369	(01) 00382902316916	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231691	Sensi Disc Piperacillin/Tazobactam – 100/10 µg	2056587	(01) 00382902316916	31-Mar-2024	single (one cartridge) or Ten Pack (10 cartridges)
231691	Sensi Disc Piperacillin/Tazobactam – 100/10 µg	2119595	(01) 00382902316916	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
231691	Sensi Disc Piperacillin/Tazobactam – 100/10 µg	3003951	(01) 00382902316916	31-Jan-2025	single (one cartridge) or Ten Pack (10 cartridges)
231691	Sensi Disc Piperacillin/Tazobactam – 100/10 µg	3062107	(01) 00382902316916	31-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
231691	Sensi Disc Piperacillin/Tazobactam – 100/10 µg	3118199	(01) 00382902316916	31-May-2025	single (one cartridge) or Ten Pack (10 cartridges)
231703	Sensi Disc Meropenem – 10 µg	2147649	(01) 0 038290 231703 6	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231705	Sensi Disc Levofloxacin – 5 µg	1334510	(01) 0 038290 231705 0	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231705	Sensi Disc Levofloxacin – 5 µg	2004644	(01) 0 038290 231705 0	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
231705	Sensi Disc Levofloxacin – 5 µg	2091272	(01) 0 038290 231705 0	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231705	Sensi Disc Levofloxacin – 5 µg	2277179	(01) 0 038290 231705 0	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
231705	Sensi Disc Levofloxacin – 5 µg	3030830	(01) 0 038290 231705 0	28-Feb-2025	single (one cartridge) or Ten Pack (10 cartridges)
231705	Sensi Disc Levofloxacin – 5 µg	3094117	(01) 0 038290 231705 0	30-Apr-2025	single (one cartridge) or Ten Pack (10 cartridges)
232174	Sensi Disc Ertapenem – 10 µg	1341479	(01) 0 038290 232174 3	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
232174	Sensi Disc Ertapenem – 10 µg	2119610	(01) 0 038290 232174 3	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
232174	Sensi Disc Ertapenem – 10 µg	2147329	(01) 0 038290 232174 3	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
232174	Sensi Disc Ertapenem – 10 µg	2242618	(01) 0 038290 232174 3	30-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
232174	Sensi Disc Ertapenem – 10 µg	2276113	(01) 0 038290 232174 3	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
232174	Sensi Disc Ertapenem – 10 µg	2333166	(01) 0 038290 232174 3	31-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
232174	Sensi Disc Ertapenem – 10 µg	3152004	(01) 0 038290 232174 3	08-Jun-2025	single (one cartridge) or Ten Pack (10 cartridges)
232174	Sensi Disc Ertapenem – 10 µg	3181140	(01) 0 038290 232174 3	29-Jun-2025	single (one cartridge) or Ten Pack (10 cartridges)
291015	Disc Augmentin Amc-30	3003894	(01) 00382902910152	31-Jul-2024	single (one cartridge) or Ten Pack (10 cartridges)
291015	Disc Augmentin Amc-30	3271840	(01) 00382902910152	29-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)



2 International Business Park Road  
 The Strategy #08-08  
 Singapore 609930  
 Registration No. 201114149N  
[bd.com](http://bd.com)

Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
291048	Disc Cefaclor Cec-30	2333108	(01) 00382902910480	31-Dec-2025	single (one cartridge) or Ten Pack (10 cartridges)
291085	Disc Ceftazidime Caz-30	3214798	(01) 00382902910855	31-Jan-2025	single (one cartridge) or Ten Pack (10 cartridges)
291085	Disc Ceftazidime Caz-30	3271784	(01) 00382902910855	29-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
291148	Disc Imipenen Ipm-10	3118317	(01) 00382902911487	31-May-2025	single (one cartridge) or Ten Pack (10 cartridges)
291169	Disc Cefixime Cfm-5	3248087	(01) 00382902911692	04-Sep-2025	single (one cartridge) or Ten Pack (10 cartridges)
291240	Disc Ampicillin/Sulbactam Sam-20 1Ea	3040399	(01) 00382902912408	31-Aug-2024	single (one cartridge) or Ten Pack (10 cartridges)
291373	Disc Clarithromycin Clr-15	2301505	(01) 00382902913733	30-Nov-2024	single (one cartridge) or Ten Pack (10 cartridges)
291435	Disc Levofloxacin Lvx-5	3118306	(01) 00382902914358	31-May-2025	single (one cartridge) or Ten Pack (10 cartridges)
291435	Disc Levofloxacin Lvx-5	3118306	(01) 00382902914358	31-May-2025	single (one cartridge) or Ten Pack (10 cartridges)
296641	Disc Sulfamethoxazole W/Trim Sxt	3214768	(01) 00382902966418	01-Aug-2026	single (one cartridge) or Ten Pack (10 cartridges)
230733	Sensi Disc Chloramphenicol - 30 µg	1020750	(01) 0 038290 230733 4	31-Jan-2025	single (one cartridge) or Ten Pack (10 cartridges)
230733	Sensi Disc Chloramphenicol - 30 µg	1116997	(01) 0 038290 230733 4	31-May-2025	single (one cartridge) or Ten Pack (10 cartridges)
230733	Sensi Disc Chloramphenicol - 30 µg	1180198	(01) 0 038290 230733 4	31-Jul-2025	single (one cartridge) or Ten Pack (10 cartridges)
230733	Sensi Disc Chloramphenicol - 30 µg	1211448	(01) 0 038290 230733 4	30-Nov-2025	single (one cartridge) or Ten Pack (10 cartridges)
230733	Sensi Disc Chloramphenicol - 30 µg	1277739	(01) 0 038290 230733 4	31-Oct-2025	single (one cartridge) or Ten Pack (10 cartridges)
230733	Sensi Disc Chloramphenicol - 30 µg	1301482	(01) 0 038290 230733 4	30-Nov-2025	single (one cartridge) or Ten Pack (10 cartridges)
230733	Sensi Disc Chloramphenicol - 30 µg	1334126	(01) 0 038290 230733 4	31-Dec-2025	single (one cartridge) or Ten Pack (10 cartridges)
230733	Sensi Disc Chloramphenicol - 30 µg	2025900	(01) 0 038290 230733 4	28-Feb-2026	single (one cartridge) or Ten Pack (10 cartridges)
230733	Sensi Disc Chloramphenicol - 30 µg	2063283	(01) 0 038290 230733 4	31-Mar-2026	single (one cartridge) or Ten Pack (10 cartridges)
230733	Sensi Disc Chloramphenicol - 30 µg	2210054	(01) 0 038290 230733 4	31-Aug-2026	single (one cartridge) or Ten Pack (10 cartridges)
230733	Sensi Disc Chloramphenicol - 30 µg	2277205	(01) 0 038290 230733 4	31-Oct-2026	single (one cartridge) or Ten Pack (10 cartridges)
230733	Sensi Disc Chloramphenicol - 30 µg	3003821	(01) 0 038290 230733 4	31-Jan-2027	single (one cartridge) or Ten Pack (10 cartridges)
230733	Sensi Disc Chloramphenicol - 30 µg	3129527	(01) 0 038290 230733 4	31-May-2027	single (one cartridge) or Ten Pack (10 cartridges)
230733	Sensi Disc Chloramphenicol - 30 µg	9336993	(01) 0 038290 230733 4	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
230733	Sensi Disc Chloramphenicol - 30 µg	0080290	(01) 0 038290 230733 4	31-Mar-2024	single (one cartridge) or Ten Pack (10 cartridges)



2 International Business Park Road  
 The Strategy #08-08  
 Singapore 609930  
 Registration No. 201114149N  
**bd.com**

Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
230733	Sensi Disc Chloramphenicol - 30 µg	0178844	(01) 0 038290 230733 4	31-Jul-2024	single (one cartridge) or Ten Pack (10 cartridges)
230733	Sensi Disc Chloramphenicol - 30 µg	0258758	(01) 0 038290 230733 4	30-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
231251	Sensi Disc Minocycline – 30 µg	1007140	(01) 3 038290 231251 3	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
231251	Sensi Disc Minocycline – 30 µg	2146396	(01) 3 038290 231251 3	30-Jun-2025	single (one cartridge) or Ten Pack (10 cartridges)
231251	Sensi Disc Minocycline – 30 µg	1085656	(01) 3 038290 231251 3	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231251	Sensi Disc Minocycline – 30 µg	1210753	(01) 3 038290 231251 3	31-Aug-2024	single (one cartridge) or Ten Pack (10 cartridges)
231251	Sensi Disc Minocycline – 30 µg	1277053	(01) 3 038290 231251 3	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
231251	Sensi Disc Minocycline – 30 µg	1301694	(01) 3 038290 231251 3	30-Nov-2024	single (one cartridge) or Ten Pack (10 cartridges)
231251	Sensi Disc Minocycline – 30 µg	1334451	(01) 3 038290 231251 3	31-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
231251	Sensi Disc Minocycline – 30 µg	2025108	(01) 3 038290 231251 3	28-Feb-2025	single (one cartridge) or Ten Pack (10 cartridges)
231251	Sensi Disc Minocycline – 30 µg	2119580	(01) 3 038290 231251 3	31-May-2025	single (one cartridge) or Ten Pack (10 cartridges)
231263	Sensi Disc Ampicillin - 2 µg	3010977	(01) 3 038290 231263 6	31-Jan-2025	single (one cartridge) or Ten Pack (10 cartridges)
231263	Sensi Disc Ampicillin - 2 µg	3058508	(01) 3 038290 231263 6	31-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
231263	Sensi Disc Ampicillin - 2 µg	3234190	(01) 3 038290 231263 6	21-Aug-2025	single (one cartridge) or Ten Pack (10 cartridges)
231263	Sensi Disc Ampicillin - 2 µg	2003517	(01) 3 038290 231263 6	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
231263	Sensi Disc Ampicillin - 2 µg	2146414	(01) 3 038290 231263 6	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
231263	Sensi Disc Ampicillin - 2 µg	2243865	(01) 3 038290 231263 6	30-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
231263	Sensi Disc Ampicillin - 2 µg	2277199	(01) 3 038290 231263 6	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
231263	Sensi Disc Ampicillin - 2 µg	2339360	(01) 3 038290 231263 6	31-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
231264	Sensi Disc Ampicillin - 10 µg	2333100	(01) 3 038290 231264 3	31-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
231264	Sensi Disc Ampicillin - 10 µg	3030225	(01) 3 038290 231264 3	28-Feb-2025	single (one cartridge) or Ten Pack (10 cartridges)
231264	Sensi Disc Ampicillin - 10 µg	1334111	(01) 3 038290 231264 3	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231264	Sensi Disc Ampicillin - 10 µg	2003260	(01) 3 038290 231264 3	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
231264	Sensi Disc Ampicillin - 10 µg	2025840	(01) 3 038290 231264 3	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231264	Sensi Disc Ampicillin - 10 µg	2146318	(01) 3 038290 231264 3	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)





2 International Business Park Road  
 The Strategy #08-08  
 Singapore 609930  
 Registration No. 201114149N  
**bd.com**

Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
231264	Sensi Disc Ampicillin - 10 µg	2242344	(01) 3 038290 231264 3	30-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
231264	Sensi Disc Ampicillin - 10 µg	2276051	(01) 3 038290 231264 3	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	2277205	(01) 3 038290 231274 2	31-Oct-2026	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	3003821	(01) 3 038290 231274 2	31-Jan-2027	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	3129527	(01) 3 038290 231274 2	31-May-2027	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	3263117	(01) 3 038290 231274 2	19-Sep-2027	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	0057245	(01) 3 038290 231274 2	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	0080290	(01) 3 038290 231274 2	31-Mar-2024	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	1057453	(01) 3 038290 231274 2	31-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	1116997	(01) 3 038290 231274 2	31-May-2025	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	1148461	(01) 3 038290 231274 2	30-Jun-2025	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	1180198	(01) 3 038290 231274 2	31-Jul-2025	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	1211448	(01) 3 038290 231274 2	31-Aug-2025	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	1277739	(01) 3 038290 231274 2	31-Oct-2025	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	1301482	(01) 3 038290 231274 2	30-Nov-2025	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	1334126	(01) 3 038290 231274 2	31-Dec-2025	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	2025900	(01) 3 038290 231274 2	28-Feb-2026	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	2063283	(01) 3 038290 231274 2	31-Mar-2026	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	2210054	(01) 3 038290 231274 2	31-Aug-2026	single (one cartridge) or Ten Pack (10 cartridges)
231274	Sensi Disc Chloramphenicol - 30 µg	9336993	(01) 3 038290 231274 2	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	3010027	(01) 3 038290 231344 2	31-Jan-2028	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	3062263	(01) 3 038290 231344 2	31-Mar-2028	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	3094187	(01) 3 038290 231344 2	30-Apr-2028	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	3158029	(01) 3 038290 231344 2	07-Jun-2028	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	3212760	(01) 3 038290 231344 2	29-Jul-2028	single (one cartridge) or Ten Pack (10 cartridges)



Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
231344	Sensi Disc Tetracycline – 30 µg	0065366	(01) 3 038290 231344 2	31-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	0085202	(01) 3 038290 231344 2	30-Apr-2025	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	0230453	(01) 3 038290 231344 2	31-Aug-2025	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	0274792	(01) 3 038290 231344 2	31-Oct-2025	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	0335906	(01) 3 038290 231344 2	31-Dec-2025	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	1095603	(01) 3 038290 231344 2	30-Apr-2026	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	1148539	(01) 3 038290 231344 2	30-Jun-2026	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	1180271	(01) 3 038290 231344 2	31-Jul-2026	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	1242880	(01) 3 038290 231344 2	30-Sep-2026	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	1277918	(01) 3 038290 231344 2	31-Oct-2026	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	1301570	(01) 3 038290 231344 2	30-Nov-2026	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	2025088	(01) 3 038290 231344 2	28-Feb-2027	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	2056233	(01) 3 038290 231344 2	31-Mar-2027	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	2089874	(01) 3 038290 231344 2	30-Apr-2027	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	2118569	(01) 3 038290 231344 2	31-May-2027	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	2276681	(01) 3 038290 231344 2	31-Oct-2027	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	8355562	(01) 3 038290 231344 2	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	9009914	(01) 3 038290 231344 2	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	9046895	(01) 3 038290 231344 2	31-Mar-2024	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	9130797	(01) 3 038290 231344 2	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	9239639	(01) 3 038290 231344 2	31-Aug-2024	single (one cartridge) or Ten Pack (10 cartridges)
231344	Sensi Disc Tetracycline – 30 µg	9263061	(01) 3 038290 231344 2	30-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	2339393	(01) 3 038290 231539 2	31-Dec-2025	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	3003869	(01) 3 038290 231539 2	31-Jan-2026	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	3062089	(01) 3 038290 231539 2	31-Mar-2026	single (one cartridge) or Ten Pack (10 cartridges)



Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	3118079	(01) 3 038290 231539 2	31-May-2026	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	3248444	(01) 3 038290 231539 2	04-Sep-2026	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	0343394	(01) 3 038290 231539 2	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	1060689	(01) 3 038290 231539 2	31-Mar-2024	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	1116935	(01) 3 038290 231539 2	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	1148543	(01) 3 038290 231539 2	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	1210905	(01) 3 038290 231539 2	31-Aug-2024	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	1242939	(01) 3 038290 231539 2	30-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	1306103	(01) 3 038290 231539 2	30-Nov-2024	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	1334461	(01) 3 038290 231539 2	31-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	2004540	(01) 3 038290 231539 2	31-Jan-2025	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	2025116	(01) 3 038290 231539 2	28-Feb-2025	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	2056456	(01) 3 038290 231539 2	31-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	2108378	(01) 3 038290 231539 2	30-Apr-2025	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	2146453	(01) 3 038290 231539 2	30-Jun-2025	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	2277321	(01) 3 038290 231539 2	31-Oct-2025	single (one cartridge) or Ten Pack (10 cartridges)
231539	Sensi Disc Sulfamethoxazole with Trimethoprim 23.75/1.25 µg	2301307	(01) 3 038290 231539 2	30-Nov-2025	single (one cartridge) or Ten Pack (10 cartridges)
231544	Sensi Disc Rifampin – 5 µg	1116055	(01) 3 038290 231544 6	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
231544	Sensi Disc Rifampin – 5 µg	1211489	(01) 3 038290 231544 6	31-Aug-2024	single (one cartridge) or Ten Pack (10 cartridges)
231544	Sensi Disc Rifampin – 5 µg	2339407	(01) 3 038290 231544 6	31-Dec-2025	single (one cartridge) or Ten Pack (10 cartridges)
231544	Sensi Disc Rifampin – 5 µg	3030278	(01) 3 038290 231544 6	28-Feb-2026	single (one cartridge) or Ten Pack (10 cartridges)
231544	Sensi Disc Rifampin – 5 µg	3158059	(01) 3 038290 231544 6	08-Jun-2026	single (one cartridge) or Ten Pack (10 cartridges)
231544	Sensi Disc Rifampin – 5 µg	0343413	(01) 3 038290 231544 6	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231544	Sensi Disc Rifampin – 5 µg	1032232	(01) 3 038290 231544 6	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231544	Sensi Disc Rifampin – 5 µg	1088790	(01) 3 038290 231544 6	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)



Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
231544	Sensi Disc Rifampin – 5 µg	1148547	(01) 3 038290 231544 6	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
231544	Sensi Disc Rifampin – 5 µg	1278110	(01) 3 038290 231544 6	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
231544	Sensi Disc Rifampin – 5 µg	1334479	(01) 3 038290 231544 6	31-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
231544	Sensi Disc Rifampin – 5 µg	2004554	(01) 3 038290 231544 6	31-Jan-2025	single (one cartridge) or Ten Pack (10 cartridges)
231544	Sensi Disc Rifampin – 5 µg	2025125	(01) 3 038290 231544 6	28-Feb-2025	single (one cartridge) or Ten Pack (10 cartridges)
231544	Sensi Disc Rifampin – 5 µg	2056466	(01) 3 038290 231544 6	31-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
231544	Sensi Disc Rifampin – 5 µg	2090443	(01) 3 038290 231544 6	30-Apr-2025	single (one cartridge) or Ten Pack (10 cartridges)
231544	Sensi Disc Rifampin – 5 µg	2242510	(01) 3 038290 231544 6	30-Sep-2025	single (one cartridge) or Ten Pack (10 cartridges)
231607	Sensi Disc Cefotaxime – 30 µg	2301360	(01) 3 038290 231607 8	30-Nov-2024	single (one cartridge) or Ten Pack (10 cartridges)
231607	Sensi Disc Cefotaxime – 30 µg	3062270	(01) 3 038290 231607 8	31-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
231607	Sensi Disc Cefotaxime – 30 µg	3158078	(01) 3 038290 231607 8	08-Jun-2025	single (one cartridge) or Ten Pack (10 cartridges)
231607	Sensi Disc Cefotaxime – 30 µg	3261093	(01) 3 038290 231607 8	17-Sep-2025	single (one cartridge) or Ten Pack (10 cartridges)
231607	Sensi Disc Cefotaxime – 30 µg	2004580	(01) 3 038290 231607 8	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
231607	Sensi Disc Cefotaxime – 30 µg	2025164	(01) 3 038290 231607 8	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231607	Sensi Disc Cefotaxime – 30 µg	2090467	(01) 3 038290 231607 8	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231607	Sensi Disc Cefotaxime – 30 µg	2146516	(01) 3 038290 231607 8	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
231607	Sensi Disc Cefotaxime – 30 µg	2180432	(01) 3 038290 231607 8	31-Jul-2024	single (one cartridge) or Ten Pack (10 cartridges)
231607	Sensi Disc Cefotaxime – 30 µg	3003886	(01) 3 038290 231607 8	31-Jan-2025	single (one cartridge) or Ten Pack (10 cartridges)
231621	Sensi Disc Cefuroxime – 30 µg	2304487	(01) 3 038290 2316214	30-Nov-2024	single (one cartridge) or Ten Pack (10 cartridges)
231621	Sensi Disc Cefuroxime – 30 µg	3093372	(01) 3 038290 2316214	30-Apr-2025	single (one cartridge) or Ten Pack (10 cartridges)
231621	Sensi Disc Cefuroxime – 30 µg	3158087	(01) 3 038290 2316214	25-Aug-2025	single (one cartridge) or Ten Pack (10 cartridges)
231621	Sensi Disc Cefuroxime – 30 µg	3223379	(01) 3 038290 2316214	10-Aug-2025	single (one cartridge) or Ten Pack (10 cartridges)
231621	Sensi Disc Cefuroxime – 30 µg	1335657	(01) 3 038290 2316214	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231621	Sensi Disc Cefuroxime – 30 µg	2004439	(01) 3 038290 2316214	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
231621	Sensi Disc Cefuroxime – 30 µg	2090487	(01) 3 038290 2316214	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)



Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
231621	Sensi Disc Cefuroxime – 30 µg	2147701	(01) 3 038290 2316214	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
231621	Sensi Disc Cefuroxime – 30 µg	2210631	(01) 3 038290 2316214	31-Aug-2024	single (one cartridge) or Ten Pack (10 cartridges)
231629	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	3163306	(01) 3 038290 231629 0	11-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
231629	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	3163307	(01) 3 038290 231629 0	11-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
231629	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	3163310	(01) 3 038290 231629 0	11-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
231629	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	3254307	(01) 3 038290 231629 0	12-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
231629	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	2146520	(01) 3 038290 231629 0	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231629	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	2146526	(01) 3 038290 231629 0	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231629	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	2147225	(01) 3 038290 231629 0	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231629	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	2210201	(01) 3 038290 231629 0	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231629	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	2224295	(01) 3 038290 231629 0	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231629	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	2276143	(01) 3 038290 231629 0	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231629	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	2276144	(01) 3 038290 231629 0	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231629	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	3003894	(01) 3 038290 231629 0	31-Jul-2024	single (one cartridge) or Ten Pack (10 cartridges)
231629	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	3062095	(01) 3 038290 231629 0	30-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
231629	Sensi Disc Amoxicillin with Clavulanic Acid – 20/10 µg	3090649	(01) 3 038290 231629 0	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
231632	Sensi Disc Ceftazidime - 30 µg	2004587	(01) 0 038290 231632 9	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
231632	Sensi Disc Ceftazidime - 30 µg	2028883	(01) 0 038290 231632 9	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231632	Sensi Disc Ceftazidime - 30 µg	2119589	(01) 0 038290 231632 9	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
231632	Sensi Disc Ceftazidime - 30 µg	2146527	(01) 0 038290 231632 9	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
231632	Sensi Disc Ceftazidime - 30 µg	2242568	(01) 0 038290 231632 9	30-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
231632	Sensi Disc Ceftazidime - 30 µg	2301371	(01) 0 038290 231632 9	30-Nov-2024	single (one cartridge) or Ten Pack (10 cartridges)
231633	Sensi Disc Ceftazidime - 30 µg	2242568	(01) 3 038290 231633 7	30-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
231633	Sensi Disc Ceftazidime - 30 µg	2301371	(01) 3 038290 231633 7	30-Nov-2024	single (one cartridge) or Ten Pack (10 cartridges)
231633	Sensi Disc Ceftazidime - 30 µg	3062286	(01) 3 038290 231633 7	31-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)



Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
231633	Sensi Disc Ceftriaxone - 30 µg	3118111	(01) 3 038290 231633 7	31-May-2025	single (one cartridge) or Ten Pack (10 cartridges)
231633	Sensi Disc Ceftriaxone - 30 µg	3181037	(01) 3 038290 231633 7	29-Jun-2025	single (one cartridge) or Ten Pack (10 cartridges)
231633	Sensi Disc Ceftriaxone - 30 µg	2004587	(01) 3 038290 231633 7	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
231633	Sensi Disc Ceftriaxone - 30 µg	2028232	(01) 3 038290 231633 7	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231633	Sensi Disc Ceftriaxone - 30 µg	2028883	(01) 3 038290 231633 7	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231633	Sensi Disc Ceftriaxone - 30 µg	2028894	(01) 3 038290 231633 7	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231633	Sensi Disc Ceftriaxone - 30 µg	2056528	(01) 3 038290 231633 7	31-Mar-2024	single (one cartridge) or Ten Pack (10 cartridges)
231633	Sensi Disc Ceftriaxone - 30 µg	2119589	(01) 3 038290 231633 7	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
231633	Sensi Disc Ceftriaxone - 30 µg	2146527	(01) 3 038290 231633 7	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
231635	Sensi Disc Ceftriaxone - 30 µg	3062288	(01) 3 038290 231635 1	31-Mar-2026	single (one cartridge) or Ten Pack (10 cartridges)
231635	Sensi Disc Ceftriaxone - 30 µg	3158099	(01) 3 038290 231635 1	08-Jun-2026	single (one cartridge) or Ten Pack (10 cartridges)
231635	Sensi Disc Ceftriaxone - 30 µg	0344489	(01) 3 038290 231635 1	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231635	Sensi Disc Ceftriaxone - 30 µg	1033870	(01) 3 038290 231635 1	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231635	Sensi Disc Ceftriaxone - 30 µg	1088039	(01) 3 038290 231635 1	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231635	Sensi Disc Ceftriaxone - 30 µg	1116107	(01) 3 038290 231635 1	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
231635	Sensi Disc Ceftriaxone - 30 µg	1148606	(01) 3 038290 231635 1	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
231635	Sensi Disc Ceftriaxone - 30 µg	1180379	(01) 3 038290 231635 1	31-Jul-2024	single (one cartridge) or Ten Pack (10 cartridges)
231635	Sensi Disc Ceftriaxone - 30 µg	1211608	(01) 3 038290 231635 1	31-Aug-2024	single (one cartridge) or Ten Pack (10 cartridges)
231635	Sensi Disc Ceftriaxone - 30 µg	1244758	(01) 3 038290 231635 1	30-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
231635	Sensi Disc Ceftriaxone - 30 µg	1335697	(01) 3 038290 231635 1	31-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
231635	Sensi Disc Ceftriaxone - 30 µg	2026381	(01) 3 038290 231635 1	28-Feb-2025	single (one cartridge) or Ten Pack (10 cartridges)
231635	Sensi Disc Ceftriaxone - 30 µg	2028922	(01) 3 038290 231635 1	28-Feb-2025	single (one cartridge) or Ten Pack (10 cartridges)
231635	Sensi Disc Ceftriaxone - 30 µg	2056573	(01) 3 038290 231635 1	31-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
231635	Sensi Disc Ceftriaxone - 30 µg	2180451	(01) 3 038290 231635 1	31-Jul-2025	single (one cartridge) or Ten Pack (10 cartridges)
231635	Sensi Disc Ceftriaxone - 30 µg	2277488	(01) 3 038290 231635 1	31-Oct-2025	single (one cartridge) or Ten Pack (10 cartridges)



Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
231635	Sensi Disc Ceftriaxone - 30 µg	2277490	(01) 3 038290 231635 1	31-Oct-2025	single (one cartridge) or Ten Pack (10 cartridges)
231640	Sensi Disc Aztreonam - 30 µg	1211613	(01) 0 038290 231640 4	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231640	Sensi Disc Aztreonam - 30 µg	1243750	(01) 0 038290 231640 4	31-Mar-2024	single (one cartridge) or Ten Pack (10 cartridges)
231640	Sensi Disc Aztreonam - 30 µg	1278287	(01) 0 038290 231640 4	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231641	Sensi Disc Aztreonam - 30 µg	2304491	(01) 3 038290 231641 2	31-May-2025	single (one cartridge) or Ten Pack (10 cartridges)
231641	Sensi Disc Aztreonam - 30 µg	3093375	(01) 3 038290 231641 2	31-Mar-2024	single (one cartridge) or Ten Pack (10 cartridges)
231641	Sensi Disc Aztreonam - 30 µg	3184116	(01) 3 038290 231641 2	01-Jan-2026	single (one cartridge) or Ten Pack (10 cartridges)
231641	Sensi Disc Aztreonam - 30 µg	3261101	(01) 3 038290 231641 2	19-Mar-2026	single (one cartridge) or Ten Pack (10 cartridges)
231641	Sensi Disc Aztreonam - 30 µg	1211613	(01) 3 038290 231641 2	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231641	Sensi Disc Aztreonam - 30 µg	1278287	(01) 3 038290 231641 2	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231641	Sensi Disc Aztreonam - 30 µg	1305148	(01) 3 038290 231641 2	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
231641	Sensi Disc Aztreonam - 30 µg	2004458	(01) 3 038290 231641 2	31-Jul-2024	single (one cartridge) or Ten Pack (10 cartridges)
231641	Sensi Disc Aztreonam - 30 µg	2028924	(01) 3 038290 231641 2	31-Aug-2024	single (one cartridge) or Ten Pack (10 cartridges)
231641	Sensi Disc Aztreonam - 30 µg	2090502	(01) 3 038290 231641 2	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
231641	Sensi Disc Aztreonam - 30 µg	2244457	(01) 3 038290 231641 2	31-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
231645	Sensi Disc Imipenem – 10 µg	3062101	(01) 3 038290 231645 0	31-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
231645	Sensi Disc Imipenem – 10 µg	3094134	(01) 3 038290 231645 0	30-Apr-2025	single (one cartridge) or Ten Pack (10 cartridges)
231645	Sensi Disc Imipenem – 10 µg	1335006	(01) 3 038290 231645 0	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231645	Sensi Disc Imipenem – 10 µg	2056580	(01) 3 038290 231645 0	31-Mar-2024	single (one cartridge) or Ten Pack (10 cartridges)
231645	Sensi Disc Imipenem – 10 µg	2081151	(01) 3 038290 231645 0	31-Mar-2024	single (one cartridge) or Ten Pack (10 cartridges)
231645	Sensi Disc Imipenem – 10 µg	2091918	(01) 3 038290 231645 0	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231645	Sensi Disc Imipenem – 10 µg	2091925	(01) 3 038290 231645 0	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231645	Sensi Disc Imipenem – 10 µg	2147704	(01) 3 038290 231645 0	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
231645	Sensi Disc Imipenem – 10 µg	2210318	(01) 3 038290 231645 0	31-Aug-2024	single (one cartridge) or Ten Pack (10 cartridges)
231645	Sensi Disc Imipenem – 10 µg	2210343	(01) 3 038290 231645 0	31-Aug-2024	single (one cartridge) or Ten Pack (10 cartridges)



Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
231645	Sensi Disc Imipenem – 10 µg	2277500	(01) 3 038290 231645 0	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
231645	Sensi Disc Imipenem – 10 µg	2301385	(01) 3 038290 231645 0	30-Nov-2024	single (one cartridge) or Ten Pack (10 cartridges)
231653	Sensi Disc Cefaclor - 30 µg	1085898	(01) 3 038290 231653 5	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231653	Sensi Disc Cefaclor - 30 µg	1273095	(01) 3 038290 231653 5	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
231653	Sensi Disc Cefaclor - 30 µg	1273096	(01) 3 038290 231653 5	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
231653	Sensi Disc Cefaclor - 30 µg	2004485	(01) 3 038290 231653 5	31-Jan-2025	single (one cartridge) or Ten Pack (10 cartridges)
231653	Sensi Disc Cefaclor - 30 µg	2147706	(01) 3 038290 231653 5	30-Jun-2025	single (one cartridge) or Ten Pack (10 cartridges)
231653	Sensi Disc Cefaclor - 30 µg	2333108	(01) 3 038290 231653 5	31-Dec-2025	single (one cartridge) or Ten Pack (10 cartridges)
231657	Sensi Disc Ciprofloxacin – 5 µg	1060926	(01) 00382902316572	31-Mar-2024	single (one cartridge) or Ten Pack (10 cartridges)
231657	Sensi Disc Ciprofloxacin – 5 µg	1279842	(01) 00382902316572	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
231657	Sensi Disc Ciprofloxacin – 5 µg	2028533	(01) 00382902316572	28-Feb-2025	single (one cartridge) or Ten Pack (10 cartridges)
231657	Sensi Disc Ciprofloxacin – 5 µg	2091948	(01) 00382902316572	30-Apr-2025	single (one cartridge) or Ten Pack (10 cartridges)
231657	Sensi Disc Ciprofloxacin – 5 µg	2301402	(01) 00382902316572	30-Nov-2025	single (one cartridge) or Ten Pack (10 cartridges)
231658	Sensi Disc Ciprofloxacin – 5 µg	1335458	(01) 30382902316580	31-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
231658	Sensi Disc Ciprofloxacin – 5 µg	2336475	(01) 30382902316580	31-Dec-2025	single (one cartridge) or Ten Pack (10 cartridges)
231658	Sensi Disc Ciprofloxacin – 5 µg	3062295	(01) 30382902316580	31-Mar-2026	single (one cartridge) or Ten Pack (10 cartridges)
231658	Sensi Disc Ciprofloxacin – 5 µg	3093427	(01) 30382902316580	30-Apr-2026	single (one cartridge) or Ten Pack (10 cartridges)
231658	Sensi Disc Ciprofloxacin – 5 µg	3163048	(01) 30382902316580	11-Jun-2026	single (one cartridge) or Ten Pack (10 cartridges)
231658	Sensi Disc Ciprofloxacin – 5 µg	0344330	(01) 30382902316580	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231658	Sensi Disc Ciprofloxacin – 5 µg	1007112	(01) 30382902316580	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
231658	Sensi Disc Ciprofloxacin – 5 µg	1060926	(01) 30382902316580	31-Mar-2024	single (one cartridge) or Ten Pack (10 cartridges)
231658	Sensi Disc Ciprofloxacin – 5 µg	1127037	(01) 30382902316580	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
231658	Sensi Disc Ciprofloxacin – 5 µg	1279842	(01) 30382902316580	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
231658	Sensi Disc Ciprofloxacin – 5 µg	1335457	(01) 30382902316580	31-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
231658	Sensi Disc Ciprofloxacin – 5 µg	2028533	(01) 30382902316580	28-Feb-2025	single (one cartridge) or Ten Pack (10 cartridges)





Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
231658	Sensi Disc Ciprofloxacin – 5 µg	2028547	(01) 30382902316580	28-Feb-2025	single (one cartridge) or Ten Pack (10 cartridges)
231658	Sensi Disc Ciprofloxacin – 5 µg	2091948	(01) 30382902316580	30-Apr-2025	single (one cartridge) or Ten Pack (10 cartridges)
231658	Sensi Disc Ciprofloxacin – 5 µg	2091980	(01) 30382902316580	30-Apr-2025	single (one cartridge) or Ten Pack (10 cartridges)
231658	Sensi Disc Ciprofloxacin – 5 µg	2301402	(01) 30382902316580	30-Nov-2025	single (one cartridge) or Ten Pack (10 cartridges)
231658	Sensi Disc Ciprofloxacin – 5 µg	2336474	(01) 30382902316580	31-Dec-2025	single (one cartridge) or Ten Pack (10 cartridges)
231660	Sensi Disc Ampicillin with Sulbactam – 10/10 µg	3010074	(01) 3 038290 231660 3	31-Jan-2025	single (one cartridge) or Ten Pack (10 cartridges)
231660	Sensi Disc Ampicillin with Sulbactam – 10/10 µg	3093387	(01) 3 038290 231660 3	30-Apr-2025	single (one cartridge) or Ten Pack (10 cartridges)
231660	Sensi Disc Ampicillin with Sulbactam – 10/10 µg	2063349	(01) 3 038290 231660 3	31-Mar-2024	single (one cartridge) or Ten Pack (10 cartridges)
231660	Sensi Disc Ampicillin with Sulbactam – 10/10 µg	2119619	(01) 3 038290 231660 3	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
231660	Sensi Disc Ampicillin with Sulbactam – 10/10 µg	2276123	(01) 3 038290 231660 3	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
231664	Sensi Disc Cefixime - 5 µg	2004492	(01) 30382902316641	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
231664	Sensi Disc Cefixime - 5 µg	3010077	(01) 30382902316641	31-Jan-2025	single (one cartridge) or Ten Pack (10 cartridges)
231664	Sensi Disc Cefixime - 5 µg	2091014	(01) 30382902316641	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231672	Sensi Disc Ofloxacin – 5 µg	2039254	(01) 3 038290 231672 6	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231672	Sensi Disc Ofloxacin – 5 µg	2180484	(01) 3 038290 231672 6	31-Jul-2024	single (one cartridge) or Ten Pack (10 cartridges)
231672	Sensi Disc Ofloxacin – 5 µg	2304476	(01) 3 038290 231672 6	30-Nov-2024	single (one cartridge) or Ten Pack (10 cartridges)
231674	Sensi Disc Cefpodoxime – 10 µg	3157106	(01) 30382902316740	08-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
231674	Sensi Disc Cefpodoxime – 10 µg	3109152	(01) 30382902316740	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231678	Sensi Disc Clarithromycin – 15 µg	2119623	(01) 30382902316788	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
231678	Sensi Disc Clarithromycin – 15 µg	2339435	(01) 30382902316788	31-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
231678	Sensi Disc Clarithromycin – 15 µg	3129541	(01) 30382902316788	31-May-2025	single (one cartridge) or Ten Pack (10 cartridges)
231682	Sensi Disc Azithromycin -15 µg	3093408	(01) 30382902316825	30-Apr-2025	single (one cartridge) or Ten Pack (10 cartridges)
231682	Sensi Disc Azithromycin -15 µg	3163065	(01) 30382902316825	11-Jun-2025	single (one cartridge) or Ten Pack (10 cartridges)
231682	Sensi Disc Azithromycin -15 µg	2004501	(01) 30382902316825	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
231682	Sensi Disc Azithromycin -15 µg	2147712	(01) 30382902316825	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)



Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
231682	Sensi Disc Azithromycin -15 µg	2304517	(01) 3 038290 2316825	30-Nov-2024	single (one cartridge) or Ten Pack (10 cartridges)
231692	Sensi Disc Piperacillin/Tazobactam – 100/10 µg	3003951	(01) 3 038290 231692 4	31-Jan-2025	single (one cartridge) or Ten Pack (10 cartridges)
231692	Sensi Disc Piperacillin/Tazobactam – 100/10 µg	3184117	(01) 3 038290 231692 4	02-Jul-2025	single (one cartridge) or Ten Pack (10 cartridges)
231692	Sensi Disc Piperacillin/Tazobactam – 100/10 µg	3223408	(01) 3 038290 231692 4	10-Aug-2025	single (one cartridge) or Ten Pack (10 cartridges)
231692	Sensi Disc Piperacillin/Tazobactam – 100/10 µg	2026369	(01) 3 038290 231692 4	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231692	Sensi Disc Piperacillin/Tazobactam – 100/10 µg	2056587	(01) 3 038290 231692 4	31-Mar-2024	single (one cartridge) or Ten Pack (10 cartridges)
231692	Sensi Disc Piperacillin/Tazobactam – 100/10 µg	2119595	(01) 3 038290 231692 4	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
231692	Sensi Disc Piperacillin/Tazobactam – 100/10 µg	2147716	(01) 3 038290 231692 4	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
231692	Sensi Disc Piperacillin/Tazobactam – 100/10 µg	2244463	(01) 3 038290 231692 4	30-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
231692	Sensi Disc Piperacillin/Tazobactam – 100/10 µg	3062107	(01) 3 038290 231692 4	31-Mar-2025	single (one cartridge) or Ten Pack (10 cartridges)
231692	Sensi Disc Piperacillin/Tazobactam – 100/10 µg	3118199	(01) 3 038290 231692 4	31-May-2025	single (one cartridge) or Ten Pack (10 cartridges)
231696	Sensi Disc Cefepime - 30 µg	2339439	(01) 3 038290 231696 2	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
231696	Sensi Disc Cefepime - 30 µg	3118297	(01) 3 038290 231696 2	30-Nov-2024	single (one cartridge) or Ten Pack (10 cartridges)
231696	Sensi Disc Cefepime - 30 µg	3163524	(01) 3 038290 231696 2	11-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
231696	Sensi Disc Cefepime - 30 µg	3282873	(01) 3 038290 231696 2	09-Apr-2025	single (one cartridge) or Ten Pack (10 cartridges)
231696	Sensi Disc Cefepime - 30 µg	2147720	(01) 3 038290 231696 2	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231696	Sensi Disc Cefepime - 30 µg	2213567	(01) 3 038290 231696 2	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231696	Sensi Disc Cefepime - 30 µg	2273732	(01) 3 038290 231696 2	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231704	Sensi Disc Meropenem – 10 µg	2213585	(01) 3 038290 231704 4	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231704	Sensi Disc Meropenem – 10 µg	3061281	(01) 3 038290 231704 4	30-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
231704	Sensi Disc Meropenem – 10 µg	3090652	(01) 3 038290 231704 4	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
231704	Sensi Disc Meropenem – 10 µg	3213825	(01) 3 038290 231704 4	30-Jan-2025	single (one cartridge) or Ten Pack (10 cartridges)
231704	Sensi Disc Meropenem – 10 µg	2187975	(01) 3 038290 231704 4	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
231704	Sensi Disc Meropenem – 10 µg	2213583	(01) 3 038290 231704 4	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
231706	Sensi Disc Levofloxacin – 5 µg	2277179	(01) 3 038290 231706 8	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)



Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
231706	Sensi Disc Levofloxacin – 5 µg	3030830	(01) 3 038290 231706 8	28-Feb-2025	single (one cartridge) or Ten Pack (10 cartridges)
231706	Sensi Disc Levofloxacin – 5 µg	3094117	(01) 3 038290 231706 8	30-Apr-2025	single (one cartridge) or Ten Pack (10 cartridges)
231706	Sensi Disc Levofloxacin – 5 µg	1334510	(01) 3 038290 231706 8	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
231706	Sensi Disc Levofloxacin – 5 µg	2004644	(01) 3 038290 231706 8	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
231706	Sensi Disc Levofloxacin – 5 µg	2091272	(01) 3 038290 231706 8	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
231706	Sensi Disc Levofloxacin – 5 µg	2213599	(01) 3 038290 231706 8	31-Aug-2024	single (one cartridge) or Ten Pack (10 cartridges)
231758	Sensi Disc Moxifloxacin – 5 µg	3003061	(01) 3 038290 231758 7	31-Jan-2025	single (one cartridge) or Ten Pack (10 cartridges)
231758	Sensi Disc Moxifloxacin – 5 µg	3243078	(01) 3 038290 231758 7	30-Aug-2025	single (one cartridge) or Ten Pack (10 cartridges)
231758	Sensi Disc Moxifloxacin – 5 µg	2119637	(01) 3 038290 231758 7	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
231758	Sensi Disc Moxifloxacin – 5 µg	2153502	(01) 3 038290 231758 7	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
231758	Sensi Disc Moxifloxacin – 5 µg	2301547	(01) 3 038290 231758 7	30-Nov-2024	single (one cartridge) or Ten Pack (10 cartridges)
232175	Sensi Disc Ertapenem – 10 µg	2276113	(01) 3 038290 232175 1	31-Oct-2024	single (one cartridge) or Ten Pack (10 cartridges)
232175	Sensi Disc Ertapenem – 10 µg	2333166	(01) 3 038290 232175 1	31-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
232175	Sensi Disc Ertapenem – 10 µg	2333175	(01) 3 038290 232175 1	31-Dec-2024	single (one cartridge) or Ten Pack (10 cartridges)
232175	Sensi Disc Ertapenem – 10 µg	1341479	(01) 3 038290 232175 1	31-Dec-2023	single (one cartridge) or Ten Pack (10 cartridges)
232175	Sensi Disc Ertapenem – 10 µg	2119610	(01) 3 038290 232175 1	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
232175	Sensi Disc Ertapenem – 10 µg	2119612	(01) 3 038290 232175 1	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
232175	Sensi Disc Ertapenem – 10 µg	2147329	(01) 3 038290 232175 1	30-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
232175	Sensi Disc Ertapenem – 10 µg	2242618	(01) 3 038290 232175 1	30-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)
232219	Sensi Disc Doripenem - 10 µg	3040749	(01) 3 038290 232219 2	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
232219	Sensi Disc Doripenem - 10 µg	3093444	(01) 3 038290 232219 2	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
232219	Sensi Disc Doripenem - 10 µg	3163515	(01) 3 038290 232219 2	11-Jun-2024	single (one cartridge) or Ten Pack (10 cartridges)
232231	Sensi Disc Ceftaroline – 30 µg	3116030	(01) 3 038290 232231 4	29-Feb-2024	single (one cartridge) or Ten Pack (10 cartridges)
232231	Sensi Disc Ceftaroline – 30 µg	3212690	(01) 3 038290 232231 4	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
291270	Sensi Disc Augmentin - 3µg	3093584	(01) 30382902912706	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)



2 International Business Park Road  
The Strategy #08-08  
Singapore 609930  
Registration No. 201114149N  
[bd.com](http://bd.com)

Catalog No	Product Name	Lot Number	UDI-DI	Expiration Date	Product Package Size
291270	Sensi Disc Augmentin - 3µg	3184136	(01) 30382902912706	02-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
291308	Sensi Disc Cefotaxime – 5 µg	3003046	(01) 30382902913086	31-Jan-2024	single (one cartridge) or Ten Pack (10 cartridges)
291308	Sensi Disc Cefotaxime – 5 µg	3093596	(01) 30382902913086	30-Apr-2024	single (one cartridge) or Ten Pack (10 cartridges)
291308	Sensi Disc Cefotaxime – 5 µg	3116037	(01) 30382902913086	31-May-2024	single (one cartridge) or Ten Pack (10 cartridges)
291308	Sensi Disc Cefotaxime – 5 µg	3248537	(01) 30382902913086	04-Sep-2024	single (one cartridge) or Ten Pack (10 cartridges)



**CUSTOMER RESPONSE FORM**  
**IDS-23-4851-FA**  
**BD BBL™ Sensi-Disc™**

Please fill in the information below so that we may acknowledge your receipt of this notification. Simply complete and return the completed form to **SEA\_Quality** [SEA\\_Quality@bd.com](mailto:SEA_Quality@bd.com) / local BD representative by **23 February 2024**.

Please tick as appropriate.

**For Removal (Recall): For Catalog #291270 only (Impacted Lot: 3093584 & 3184136)**

- I have read and understood the attached notice & will share this Product Recall with all users within my facility. I shall clarify with the appointed distributor/ BD representative for clarification (s) to this Product Recall.
- We do not have affected product(s) in our inventory.
- We have affected product(s) in inventory and will return to my distributor/local BD representative for destruction.

Product Description	Catalog No.	Lot No.	Quantity (eaches) received)	Quantity (eaches) sold	Remaining Quantity (eaches) in inventory to be *returned



**For Product Advisory (No disposition of product): All other Catalogs listed in Appendix 1**

I have read and understood the attached notice & will share this Advisory Notice with all users within my facility.

**Completed by:**

<b>Name:</b>	
<b>Signature:</b>	
<b>Date:</b>	
<b>Facility / Address / Telephone Number:</b>	